



MANIPULATIVE PERFORMANCE TEST

Procedures Handbook 4300

MAKING A DOUBLE DONUT
HOSE ROLL – TWO
PERSON

TOPIC/EQUIPMENT: Making a Double Donut Hose Roll - Two Person Method

CATEGORY: Level II

POINTS POSSIBLE: 100

TIME ALLOWED: Two Minutes

BEHAVIORAL OBJECTIVE:

Condition: One length of hose with male and female couplings

Behavior: The student will roll hose using a double donut roll, two person method

Standard: Within a minimum of 80% accuracy within two (2) minutes

MATERIALS NEEDED:

- Two (2) students
- One (1) section of hose with male and female couplings
- One (1) rubber band or other securing device
- One (1) manipulative performance test cover page
- Individual score sheets
- One (1) stop watch
- One (1) clip board
- One (1) red pen

PROCEDURES:

SCORING: Points will be deducted for each step omitted, performed improperly, or performed out of sequence. Steps designated by Pass/No Pass must be performed properly or the student fails the entire examination. If the student goes over the allotted time he/she will fail the entire examination.

SPECIAL NOTES:

4307.9



MANIPULATIVE PERFORMANCE TEST

Procedures Handbook 4300

MAKING A DOUBLE DONUT
HOSE ROLL – TWO
PERSON

DATE ___/___/___ TEST # _____ RETEST # _____ UNIT # _____

STUDENT'S NAME _____

EVALUATOR'S NAME _____

POINTS

TIME START

1. Lay hose fully extended and straight on the ground removing all twists and kinks. 10
2. Grasp either coupling and carry it towards the other coupling until couplings are even and adjacent to one another 15
3. Cross over hose every twenty (20) feet to keep rolls together. 15
4. First person rolls hose by rolling hose from fold back on itself to the couplings keeping the outer edges even. 15
5. Second person straddles hose while bent at waist and maintains five (5) foot interval while moving backwards keeping hose sections together, removing slack and anchoring cross overs to make tight rolls. 15
6. Once hose is fully rolled, lay the roll on its side and flatten the sides of the rolls 10
7. Secure rolled hose with a rubber band or other restraining device. 10
8. Store hose in designated area or compartment. 10

TIME STOP

POINTS DEDUCTED: _____:

POINTS POSSIBLE: 100

POINTS DEDUCTED: _____

4307.9



MANIPULATIVE PERFORMANCE TEST

Procedures Handbook 4300

MAKING A DOUBLE DONUT
HOSE ROLL – TWO
PERSON

FINAL SCORE:

EVALUATOR'S SIGNATURE: _____

DATE: _____

COMMENTS:
